
UM–SJTU JOINT INSTITUTE
PHYSICS LABORATORY
(VP241)

LABORATORY REPORT
EXERCISE 4
POLARIZATION OF LIGHT

CAO ZHIYUAN

NOVEMBER 1, 2019

PARTNERS:

NAME: CAO ZHIYUAN ID: 518370910030 GROUP: 19

NAME: GU ZHENG ID: 518370910190 GROUP: 19

DATE PERFORMED:

OCTOBER 25, 2019

Contents

1	OBJECTIVE	1
2	THEORETICAL BACKGROUND	1
2.1	POLARIZATION OF LIGHT	1
2.2	POLARIZER	2
2.3	MALUS' LAW	2
2.4	GENERATION OF ELLIPTICALLY AND CIRCULARLY POLARIZED LIGHT, HALF-WAVE AND QUARTER-WAVE PLATES	3
3	APPARATUS AND EXPERIMENTAL SETUP	5
4	PROCEDURES [1]	6
4.1	APPARATUS ADJUSTMENT	6
4.2	DEMONSTRATION OF MALUS' LAW	6
4.3	LINEARLY POLARIZED LIGHT AND THE HALF-WAVE PLATE	7
4.4	CIRCULARLY AND ELLIPTICALLY POLARIZED LIGHT AND THE 1/4-WAVE PLATE .	8
5	RESULTS	8
5.1	APPARATUS ADJUSTMENT	8
5.2	DEMONSTRATION OF MALUS' LAW	8
5.3	LINEARLY POLARIZED LIGHT AND THE HALF-WAVE PLATE	10
5.4	CIRCULARLY AND ELLIPTICALLY POLARIZED LIGHT AND THE 1/4-WAVE PLATE .	12
5.4.1	ROTATION ANGLE 0°	12
5.4.2	ROTATION ANGLE 20°	14
5.4.3	ROTATION ANGLE 45°	16
5.4.4	ROTATION ANGLE 70°	18
6	CONCLUSION AND DISCUSSION	19
6.1	ANSWER TO THE QUESTIONS IN LAB MANUAL	19
6.2	DISCUSSION	19
6.3	CONCLUSION	20
A	MEASUREMENT UNCERTAINTY ANALYSIS	20
A.1	DEMONSTRATION OF MALUS' LAW	20
A.2	LINEARLY POLARIZED LIGHT AND THE HALF-WAVE PLATE	22
A.3	CIRCULARLY AND ELLIPTICALLY POLARIZED LIGHT AND THE 1/4-WAVE PLATE .	22
A.3.1	ROTATION ANGLE 0°	23
A.3.2	ROTATION ANGLE 20°	24
A.3.3	ROTATION ANGLE 45°	25
A.3.4	ROTATION ANGLE 70°	25
B	DATA SHEET	26

1 OBJECTIVE

- Understand the basic properties of light, mainly the polarization phenomenon of light.
- Learn and verify Malus' law.
- Understand the roles half-wave and quarter-wave plates play in the optical systems.
- Investigate on how to generate and detect the circularly and elliptically polarized light.

2 THEORETICAL BACKGROUND

In our nature, waves exist everywhere. Generally speaking, waves can be classified into longitudinal and transverse waves based on the direction of oscillations. As for light, which can be described in terms of electromagnetic waves, since its plane of oscillations of the electric field vector is perpendicular to the direction of light propagation, light can be regarded as a transverse wave.

One of the most important characteristics of a transverse wave is the polarization phenomenon. In this lab, we focus on the polarization of light. The commonest light in our life is *natural light*, which is also called *unpolarized light*. For natural light, the incident light is a random mixture of waves with the electric field vector oscillating in all possible transverse directions. Also, the distribution of the electric field vector is uniform in all directions. On the other hand, for *polarized light* the distribution of directions is not uniform. Starting from these basic theories, research on wave polarization has a numerical application in the real world, ranging from optical measurement devices, to the realm of crystal structure.

2.1 POLARIZATION OF LIGHT

The electric field vector \vec{E} corresponds to the visible part of the spectrum in the context of electromagnetic waves and is sometimes called the light vector. It describes a propagating electric field that changes with time. In the plane perpendicular to the propagation direction of the wave, the direction of the light vector may be different. The light whose light vector maintains a certain direction of oscillation is called *linear polarization*, and the axis that defines the direction is called the polarization axis. (Fig.1)




Figure 1: (a) Linearly polarized light with the polarization axis perpendicular to the page plane. (b) Linearly polarized light with the polarization axis parallel to the page plane. [1]

There are other types of polarized light. For *circularly polarized light*, the direction of the light vector and the direction of the light propagation rotate so that the end point moves in a circular motion. For *elliptically polarized light*, the end point of direction moves along an ellipse. The Figure is shown in Fig.2.




Figure 2: Elliptically polarized light propagating in the z direction. The light is polarized in the xy plane. [1]

Although unpolarized light is emitted from ordinary light sources, it can be regarded as a mixture of linearly polarized light whose distribution of direction is statistically equal to each other. For partially polarized light, it can be seen as combining a natural light and a polarized light together, because a partial polarized light has a non-uniform distribution of directions of light vector, indicating it can be decomposed into a natural light and a polarized one.

2.2 POLARIZER

The device we commonly used to produce polarized light is called polarizer, which can also be called polaroid. According to the principle of dichroism, a selective absorption mechanism tends to allow the light polarized in a certain direction (direction of the crystal alignment) to pass through the material, while the light polarized in all other directions is absorbed [1]. A polarizer is able to transform an incident light which is originally natural into linearly polarized light.

Not only does a polarization device transform natural light into polarized light (the function of a polarizer), but it also act as a device to analyze different polarized light. In this case, polarizer acts as an analyzer.

2.3 MALUS' LAW

After the light passes through a polarization device, one noticeable effect is that the intensity of light brightness changes, and the light may become darker.

Assume that two polarization device are places together in a line (See Fig.3), where the left one acts as a polarizer and the right one acts as an analyzer, and given that the angle between their polarization axes is θ , the intensity of the light leaving the anlayzer is:

$$I_{light} = I_{light,0} \cos^2 \theta \quad (1)$$




Figure 3: Change in the brightness of the light depends on the mutual orientation of the polarizer and the analyzer. [1]

where I_{light} is the intensity of light passing through the polarizer, and $I_{light,0}$ is the intensity of light passing through the analyzer. Equation (1), discovered in 1809, is known as the Malus' law, which is named after Etienne-Louis Malus.

From the equation, we can clearly observe that if there is a single polarizer, the transmitted light intensity of an incident polarized light will change periodically when the polarizer is rotated, with a period $T = 2\pi$. For different type of polarized light, the phenomenon will also differ from each other:

- For an initially partially or elliptically polarized incident light, the minimum intensity will not reaches zero because of the geometric shape of light vector distribution.
- For a natural or circularly polarized light, the intensity will not change at all.

Based on above properties, by using a polarizer, we can tell linearly polarized light from the natural and circularly polarized light ones.

2.4 GENERATION OF ELLIPTICALLY AND CIRCULARLY POLARIZED LIGHT. HALF-WAVE AND QUARTER-WAVE PLATES

Assume that a linearly polarized light passes thorough a crystal plate whose surface is generally parallel to the optical axis, and the angle between the polarization axis and the optical axis of the plate is α . Then the linearly polarized light is decomposed into two different waves: one is e-wave whose oscillation direction is parallel to the optical axis of the plate while another is o-wave whose oscillation direction is perpendicular to the optical axis. They are named as extraordinary axis and ordinary axis respectively. In fact, the two waves propagate in the same direction, however, with a different speed. We can calculate the resulting optical path difference of a plate with thickness d as follows:

$$\Delta = (n_e - n_o)d \quad (2)$$

Hence, we are able to find out the phase difference accordingly:

$$\delta = \frac{2\pi}{\lambda}(n_e - n_o)d \quad (3)$$




Figure 4: Linearly polarized light passing through a waveplate. [1]

where δ is the phase difference, λ stands for the wavelength, n_e refers to the refractive index for the extraordinary axis, and n_o is the refractive index for the ordinary axis. A crystal with $\delta > 0$ is called positive crystal, while a $\delta < 0$ one is called negative crystal.

In Fig.4, we can see that there are two components of light vectors, which can be expressed as follows:

$$\begin{aligned} E_x &= A \sin \alpha \cos \omega t = A_o \cos \omega t \\ E_y &= A \cos \alpha \cos(\omega t + \delta) = A_e \cos(\omega t + \delta) \end{aligned}$$

To obtain the trajectory of the light vector, we eliminate the time t in the above equations, and we have:

$$\frac{E_x^2}{A_o^2} + \frac{E_y^2}{A_e^2} - 2 \frac{E_x E_y}{A_o A_e} \cos \delta = \sin^2 \delta \quad (4)$$

Note that Eq.(4) is the only valid for $\delta = \pm \pi/2$.

For some particular chosen thickness of plate, the optical path will change at the same time, and some interesting cases occur, which is discussed below:

- If $\Delta = k\lambda$ ($k = 0, 1, 2, \dots$), the phase difference $\delta = 0$. Eq.(4) can be rewritten into a simpler form:

$$E_y = \frac{A_e}{A_o} E_x$$

which is actually a linear equation. So the transmitted light is linearly polarized and oscillates in the same direction. A wave plate satisfying this condition is called a *full wave plate*. Light passing through a full-wave plate will not change its polarization state.

- If $\Delta = (2k + 1)\lambda/2$ ($k = 0, 1, 2, \dots$), the phase difference $\delta = \pi$. Eq.(4) can be rewritten into a simpler form:

$$E_y = -\frac{A_e}{A_o} E_x$$

Transmitted light is also linearly polarized with the polarization axis rotating 2α . A wave plate satisfying this condition is called a *1/2 wave plate* or a *half wave plate*. When the polarized light passes through the half-wave plate, the rotation angle of its polarization axis is 2α . Specifically, if $\alpha = \pi/4$, the transmitted light's polarization axis is vertical to that of the original light.

- If $\Delta = (2k + 1)\lambda/4$ ($k = 0, 1, 2, \dots$), the phase difference $\delta = \pm\pi/2$. Eq.(4) can be rewritten into a simpler form:

$$\frac{E_x^2}{A_o^2} + \frac{E_y^2}{A_e^2} = 1$$

As a result, the light is elliptically polarized. In this case, the very waveplate is called a *1/4-wave plate* or a *quarter-waveplate*, and is very important in optical experiments.

Specially, if $A_e = A_o = A$, then we obtain $E_x^2 + E_y^2 = A^2$. As a result, the transmitted light is circularly polarized. The polarization state after the light passes through the wave plate will vary, depending on the different angle. The scenarios are listed in Table 1:

α	state of the transmitted light	Polarization axis and the optical axis
0	linearly polarized	parallel
$\pi/2$	linearly polarized	perpendicular
$\pi/4$	circularly polarized	/
otherwise	elliptically polarize	/

Table 1: Relationship between angle α and polarization state.

3 APPARATUS AND EXPERIMENTAL SETUP

In this lab, we will use the following apparatus, which are placed on an optical bench:

- a semiconductor laser
- a tungsten iodine lamp
- a silicon photo-cell
- a UT51 digital universal meter
- two polarizers
- 1/2-wave and 1/4-wave plates
- a lens with a glass sheet

Then, the range or precision of the apparatus are listed in Table 2 below:

	Range/Precision
Angle	$\pm 2^\circ$
Current	$\pm 0.001 \mu A$

Table 2: Range and precision of the apparatus.

4 PROCEDURES [1]

4.1 APPARATUS ADJUSTMENT

- (1) First, we adjust the photo-cell by choosing the appropriate aperture. In this experiment, we use only 6.0 aperture, because it preserves the incident light intensity (shown in Fig.5). Otherwise, the intensity of light may get reduced, resulting in a zero reading on the universal meter.




Figure 5: Photo-cell. [1]

- (2) Second, with the laser fixed at one of the ends of the bench, we place the lens and the glass sheet in front of it. Then we make sure that the light passes through the center of the lens.
- (3) Third, the distance between the lens and the laser is adjusted to the focal length of the lens.
- (4) Forth, we move the glass sheet along the bench. If the size of the light spot on the glass varies significantly, repeat Step 2.
- (5) Eventually, we remove the glass sheet, and set the digital universal meter in the appropriate mode and range.

4.2 DEMONSTRATION OF MALUS' LAW

In this section, since the light source we use has already been linearly polarized source. The method we use to verify Malus' Law has some difference than the normal one. Details are shown as follows:

- (1) First, we fix the analyzer and rotate the polarizer. We will discover that the intensity of light changes periodically. Then, we need to choose the position of the polarizer where the light intensity is maximum. In this way, the output linearly polarized light is "preserved" and has the original intensity.
- (2) Second, we rotate the analyzer for 360° and observe a change in the light intensity to find the maximum electric current I_0 . Then, we remain the position where the light intensity is maximum.
- (3) Third, the angle of analyzer is 90° . At this point, the polarizing axes of the polarizer and the analyzer are perpendicular to each other.
- (4) Finally, we rotate the analyzer from 90° to 0° and record the magnitude of the current I every 5° . Record the values in a table and plot the graph I/I_0 vs. $\cos 2\theta$. After that, we perform linear fitting and compare the data with the theoretical result.

4.3 LINEARLY POLARIZED LIGHT AND THE HALF-WAVE PLATE

- (1) We set up the equipment on the optical bench as shown in Fig.6. A is the analyzer and P is the polarizer. Set the polarizing axes of A and P perpendicular to each other before placing the $1/2$ -wave plate in the apparatus; extinction of the light can be observed on screen.
- (2) After inserting the $1/2$ -wave plate, we rotate it to make the light extinction appear again and set this position as the initial position.
- (3) We rotate the $1/2$ -wave plate for $\alpha = 10^\circ$ from the initial position and the light extinction will be broken. Then rotate A to make the light extinction appear again, record the angle of rotation $\Delta\theta$ in a table.
- (4) Finally, we rotate the $1/2$ -wave plate for 10° from the previous position (now $\alpha = 10^\circ$) and repeat Step 3. Repeat this step (increase α) for 8 times. Plot the graph $\Delta\theta$ vs. θ .




Figure 6: Experimental setup for the $1/2$ -wave plate. [1]

4.4 CIRCULARLY AND ELLIPTICALLY POLARIZED LIGHT AND THE 1/4-WAVE PLATE

- (1) First, we set up the equipment on the optical bench as shown in Fig.6, and set the polarizing axes of A and P perpendicular to each other before placing the 1/4-wave plate in the apparatus; extinction of the light can be observed on screen. At this point the angle $\theta = 90^\circ$.
- (2) Second, after inserting the 1/4-wave plate, we rotate it to make the light extinction appear again and set this position as the initial position. At this point $\alpha = 0^\circ$. Rotate the 1/4-wave plate and observe the change in the light intensity.
- (3) Third, we rotate the analyzer for 360° , and record the light intensity (which is indicated by the current I) for every 360° . Record the data in a table.
- (4) We rotate the 1/4-wave plate for 20° , repeat Step 3.
- (5) We rotate the 1/4-wave plate for 45° , repeat Step 3.
- (6) After that, we rotate the 1/4-wave plate for 70° . Then rotate the analyzer and record its position and the magnitude of the current when the light intensity reaches a maximum.
- (7) Then, we plot the relation between the rotation angle of the analyzer and the light amplitude in polar coordinates. Normalize the amplitude by its maximum value. Mark the position recorded in Step 6 and compare it with the data recorded in Step 4.
- (8) Finally, we compare the result of Step 5 with that for the circular polarization. Plot a linear fit to the data when the angle is 45° .

5 RESULTS

5.1 APPARATUS ADJUSTMENT

According to what we have discussed in 4.1, we successfully adjust and warm up the apparatus we will use in this experiment.

5.2 DEMONSTRATION OF MALUS' LAW

In this section, we verify Malus' law, and the detailed results are shown in Table 3.

Uncertainty of θ is 2°

Maximum Electric Current I_0		$1.211 \pm 0.001[\mu A]$	
θ	$I[\mu A] \pm 0.001[\mu A]$	θ	$I[\mu A] \pm 0.001[\mu A]$
0°	1.211	50°	0.506
5°	1.203	55°	0.423
10°	1.178	60°	0.320
15°	1.122	65°	0.223
20°	1.123	70°	0.154
25°	1.056	75°	0.089
30°	0.951	80°	0.041
35°	0.855	85°	0.013
40°	0.736	90°	0.000
45°	0.620		

Table 3: Measurement data of Malus' law demonstration.

θ	I/I_0	u_{I/I_0}	$\cos^2\theta$	$u_{\cos^2\theta}$	θ	I/I_0	u_{I/I_0}	$\cos^2\theta$	$u_{\cos^2\theta}$
0°	1.0000	0.0012	1	0	50°	0.4178	0.0009	0.41	0.03
5°	0.9934	0.0012	0.992	0.006	55°	0.3493	0.0009	0.33	0.03
10°	0.9727	0.0012	0.970	0.012	60°	0.2642	0.0009	0.25	0.03
15°	0.9265	0.0011	0.933	0.017	65°	0.1841	0.0008	0.18	0.03
20°	0.9273	0.0011	0.88	0.02	70°	0.1272	0.0008	0.12	0.02
25°	0.8720	0.0011	0.82	0.03	75°	0.0735	0.0008	0.067	0.017
30°	0.7853	0.0010	0.75	0.03	80°	0.0339	0.0008	0.030	0.012
35°	0.7060	0.0010	0.67	0.03	85°	0.0107	0.0008	0.008	0.006
40°	0.6078	0.0009	0.59	0.03	90°	0.0000	0.0008	0	0
45°	0.5120	0.0009	0.50	0.03					

Table 4: Uncertainty of data of Malus' law demonstration.

After calculations, we obtain the data in Table 4. For sample calculation, we choose the case when $\theta = 5^\circ$, where $I_0 = 1.211 \pm 0.001[\mu A]$ and $I = 1.203 \pm 0.001[\mu A]$. We will obtain:

$$\frac{I}{I_0} = \frac{1.203}{1.211} = 0.9934 \pm 0.0012, \quad r_{u \frac{I}{I_0}} = 0.12\%$$

and

$$\cos^2\theta = \cos^2 5^\circ = 0.992 \pm 0.006, \quad r_{u \cos^2\theta} = 0.6\%$$

Hence our relative error is small and acceptable. Detailed calculations regarding uncertainties are shown in Appendix A.




Figure 7: Linear fit for I/I_0 vs. $\cos^2\theta$.

Then, with the help of Origin, we perform a linear fit, which is shown in Fig.7.

From Fig.7, we can see that the value of Pearson's R equals to 0.99929, which is close to 1, demonstrating that I/I_0 and $\cos^2\theta$ do have a linear relationship. Besides, since the slope is 1.01778, which is also close to 1, we obtain:

$$I/I_0 = \cos^2\theta$$

i.e., we have

$$I = I_0 \cos^2\theta$$

which matches the Malus' Law, with a relative uncertainty:

$$u_p = \frac{1.01778 - 1}{1} = 1.8\%$$

5.3 LINEARLY POLARIZED LIGHT AND THE HALF-WAVE PLATE

In this section, the data we recorded are shown in Table 4.

Then, we obtain the data with uncertainties, which is shown in Table 6.

Rotation angle of the 1/2-wave plate	Rotation angle of the analyzer [$^{\circ}$] ± 2 [$^{\circ}$]
initial	342
10 $^{\circ}$	327
20 $^{\circ}$	296
30 $^{\circ}$	276
40 $^{\circ}$	259
50 $^{\circ}$	237
60 $^{\circ}$	217
70 $^{\circ}$	198
80 $^{\circ}$	178
90 $^{\circ}$	158

Table 5: Measurement data of Linearly Polarized Light and the Half-wave Plate.

Now, we use Origin to perform a linear fit. The result is shown in Fig.8.




Figure 8: Linear fit for rotation angle of the analyzer vs. 1/2-wave plate.

From Fig.8, we can see that the value of Pearson's R equals to 0.99912, which is close to 1, demonstrating that θ_1 and θ_2 do have a linear relationship. Besides, since the slope is 2.06714, which is close to 2. Hence, we conclude that when a polarized light passes through a 1/2-wave plate, its polarization axis is rotated by 2α , with a relative uncertainty:

$$u_p = \frac{2.06714 - 2}{2} = 3\%$$

Rotation angle of the 1/2-wave plate	u_{θ_1}	Rotation angle of the analyzer [$^\circ$] ± 2 [$^\circ$]	u_{θ_2}
initial	2°	0	2°
10°	2°	15	3°
20°	2°	46	3°
30°	2°	66	3°
40°	2°	83	3°
50°	2°	105	3°
60°	2°	125	3°
70°	2°	144	3°
80°	2°	164	3°
90°	2°	184	3°

Table 6: Uncertainty of data of Linearly Polarized Light and the Half-wave Plate.

5.4 CIRCULARLY AND ELLIPTICALLY POLARIZED LIGHT AND THE 1/4-WAVE PLATE

5.4.1 ROTATION ANGLE 0°

In this section, the data is recorded in Table 7:

In order to calculate $\sqrt{\frac{I}{I_0}}$, we take the case when $\theta = 10^\circ$ as a sample calculation, and we will have:

$$\sqrt{\frac{I}{I_0}} = \sqrt{\frac{0.047}{1.452}} = 0.180 \pm 0.002$$

with relative uncertainty $r = 1.1\%$, which is acceptable. The detailed step for uncertainty calculations are shown in A.3.1.

Rotation angle of 1/4-wave plate: 0°					
Maximum Electric Current I_0 $1.452 \pm 0.001 \mu A$					
θ	$I[\mu A] \pm 0.001[\mu A]$	$\sqrt{\frac{I}{I_0}}$	θ	$I[\mu A] \pm 0.001[\mu A]$	$\sqrt{\frac{I}{I_0}}$
0°	0.000	0	180°	0.002	0.037
10°	0.047	0.180	190°	0.051	0.1874
20°	0.165	0.3371	200°	0.173	0.3452
30°	0.364	0.5007	210°	0.387	0.5163
40°	0.623	0.6550	220°	0.566	0.6243
50°	0.808	0.7460	230°	0.782	0.7339
60°	1.126	0.8806	240°	1.105	0.8724
70°	1.303	0.9473	250°	1.307	0.9488
80°	1.419	0.9886	260°	1.415	0.9872
90°	1.452	1.0000	270°	1.451	0.9997
100°	1.383	0.9760	280°	1.372	0.9721
110°	1.234	0.9219	290°	1.221	0.9170
120°	0.946	0.8072	300°	0.951	0.8093
130°	0.790	0.7376	310°	0.796	0.7404
140°	0.512	0.5938	320°	0.507	0.5909
150°	0.313	0.4643	330°	0.310	0.4621
160°	0.138	0.3083	340°	0.135	0.3049
170°	0.032	0.1485	350°	0.040	0.166

Table 7: Measurement data for the 1/4-wave plate (Rotation angle 0°).

Then, we are able to perform a linear fit in polar coordinates with the help of Origin. The result is shown in Fig.9. It can be observed that the transmitted light is linearly polarized.




Figure 9: Linear fit for rotation angle 0° .

5.4.2 ROTATION ANGLE 20°

In this section, the data is recorded in Table 8:

In order to calculate $\sqrt{\frac{I}{I_0}}$, we take the case when $\theta = 10^\circ$ as a sample calculation, and we will have:

$$\sqrt{\frac{I}{I_0}} = \sqrt{\frac{0.203}{1.201}} = 0.4111 \pm 0.0010$$

with relative uncertainty $r = 0.2\%$, which is acceptable. The detailed step for uncertainty calculations are shown in A.3.2.

Rotation angle of 1/4-wave plate: 20°					
Maximum Electric Current I_0 1.201 ± 0.001 μA					
θ	$I[\mu A] \pm 0.001[\mu A]$	$\sqrt{\frac{I}{I_0}}$	θ	$I[\mu A] \pm 0.001[\mu A]$	$\sqrt{\frac{I}{I_0}}$
0°	0.156	0.3604	180°	0.165	0.3707
10°	0.203	0.4111	190°	0.234	0.4414
20°	0.342	0.5336	200°	0.359	0.5467
30°	0.441	0.6060	210°	0.483	0.6342
40°	0.600	0.7068	220°	0.641	0.7306
50°	0.757	0.7939	230°	0.792	0.8121
60°	0.899	0.8652	240°	0.923	0.8767
70°	1.013	0.9184	250°	1.051	0.9355
80°	1.097	0.9557	260°	1.140	0.9743
90°	1.176	0.9895	270°	1.201	1.0000
100°	1.098	0.9562	280°	1.152	0.9794
110°	0.923	0.8767	290°	0.971	0.8992
120°	0.765	0.7981	300°	0.821	0.8268
130°	0.591	0.7015	310°	0.620	0.7185
140°	0.422	0.5928	320°	0.462	0.6202
150°	0.278	0.4811	330°	0.289	0.4905
160°	0.193	0.4009	340°	0.231	0.4386
170°	0.148	0.3510	350°	0.162	0.3673

Table 8: Measurement data for the 1/4-wave plate (Rotation angle 20°).

Then, we are able to perform a linear fit in polar coordinates with the help of Origin. The result is shown in Fig.10. It can be observed that the transmitted light is elliptically polarized.




Figure 10: Linear fit for rotation angle 20° .

5.4.3 ROTATION ANGLE 45°

In this section, the data is recorded in Table 9:

In order to calculate $\sqrt{\frac{I}{I_0}}$, we take the case when $\theta = 10^\circ$ as a sample calculation, and we will have:

$$\sqrt{\frac{I}{I_0}} = \sqrt{\frac{0.629}{0.705}} = 0.9446 \pm 0.0010$$

with relative uncertainty $r = 0.11\%$, which is acceptable. The detailed step for uncertainty calculations are shown in A.3.3.

Rotation angle of 1/4-wave plate: 45°					
Maximum Electric Current I_0 $0.705 \pm 0.001 \mu A$					
θ	$I[\mu A] \pm 0.001[\mu A]$	$\sqrt{\frac{I}{I_0}}$	θ	$I[\mu A] \pm 0.001[\mu A]$	$\sqrt{\frac{I}{I_0}}$
0°	0.612	0.9317	180°	0.9438	0.3707
10°	0.629	0.9446	190°	0.9400	0.4414
20°	0.623	0.9400	200°	0.9416	0.5467
30°	0.634	0.9483	210°	0.9506	0.6342
40°	0.650	0.9602	220°	0.9617	0.7306
50°	0.673	0.9770	230°	0.9785	0.8121
60°	0.692	0.9907	240°	0.9893	0.8767
70°	0.681	0.9828	250°	1.0000	0.9355
80°	0.692	0.9907	260°	0.9957	0.9743
90°	0.686	0.9864	270°	0.9950	1.0000
100°	0.685	0.9857	280°	0.9843	0.9794
110°	0.680	0.9821	290°	0.9763	0.8992
120°	0.675	0.9885	300°	0.9734	0.8268
130°	0.661	0.9683	310°	0.9668	0.7185
140°	0.649	0.9595	320°	0.9572	0.6202
150°	0.628	0.9438	330°	0.9438	0.4905
160°	0.620	0.9378	340°	0.9355	0.4386
170°	0.615	0.9340	350°	0.9378	0.3673

Table 9: Measurement data for the 1/4-wave plate (Rotation angle 45°).

Then, we are able to perform a linear fit in polar coordinates with the help of Origin. The result is shown in Fig.11. It can be observed that the transmitted light is circularly polarized.

Also, we apply a linear fit to the plot. We find that the value of Pearson's R is 0.02402, and Adj. R-Square is -0.02882. However, the curve is not a perfect circle. We will discuss this phenomenon later.




Figure 11: Linear fit for rotation angle 0° .

5.4.4 ROTATION ANGLE 70°

In this part, we rotate the analyzer, and find the position where the light intensity reaches a maximum. The data has been recorded in Table 10.

Rotation angle of 1/4-wave plate: 70°	
$\theta[^\circ] \pm 2[^\circ]$	268
$I[\mu A] \pm 0.001[\mu A]$	1.221

Table 10: Measurement data for the 1/4-wave plate (Rotation angle 45°).

In section 5.4.2, we have marked the point as "F" in Fig.10 in red.

6 CONCLUSION AND DISCUSSION

6.1 ANSWER TO THE QUESTIONS IN LAB MANUAL

- When the 1/2-wave plate rotates for 360° , **4** times of light extinction can be observed.
- When the analyzer rotates for 360° , **2** times of light extinction can be observed.
- The polarization state won't change. Namely, the transmitted light is still linearly polarized. This is because the function of the 1/2-wave plate is simply to make the polarization axis rotated by angle 2α . For example, if $\alpha = 45^\circ$, then the polarization axis of the transmitted light is vertical to that of the original light.

6.2 DISCUSSION

In this lab, generally speaking, we do a great job and have some satisfactory results. We verify Malus' Law, and focus on the properties of 1/2 and 1/4-wave plates. The uncertainty is not so large. However, there still exists some errors, like the polar linear fit has a smaller Pearson's r . Hence, some possible sources that may lead to errors are listed below:

- Although when doing this lab, we are in a dark environment, light from other groups may disturb us and influence our results.
- When doing the lab, readings on the digital meter "oscillate" for a long time. In other words, it's hard for the digital meter to be in a steady to ensure an accurate value.
- The polarizer and analyzer may not be parallel in a line with each other, leading to deviation in refraction.
- The polarizer's neck is not tight enough, which can be rotated for some small angle.
- Sometimes in this lab, because of darkness, we may misread some value, or rotate the polarizer unintentionally.

Based on the above factors and my own thoughts, some suggestions to improve the results are listed below:

- We can add some window shade between each group so that the light emitted from one group will not influence another group.
- The design of the polarizer can be modified. To be more specific, a buckle can be added on the polarizer to prevent it from rotating when we don't rotate it intentionally.
- For oscillating values on the digital meter, we can choose the value in the middle, or we can wait for a longer time for it to be steady. Also, we can do multiple times of experiment and take their average to minimize the errors.
- The distance between the polarizer and analyzer should be appropriate. Thus we avoid to have a too large or too small current.

- For simplicity and accuracy, when rotating the polarizer and analyzer, we should always rotate it in one specific direction.
- For some old devices, we can change them into new ones.

6.3 CONCLUSION

In this lab, we learned a lot, including:

- We understand the polarization phenomenon of light.
- We study and verify Malus' Law.
- We get familiar with 1/2 and 1/4-wave plate which is important in optical systems.
- Generate, analyze and identify the linearly polarized, circularly polarized, and elliptically polarized light.
- Learn how to use UT51 digital universal meter, semiconductor laser, and various useful optical apparatuses.

To sum up, in this lab we really learned a lot. Not only do we improve our cooperating skills, but we analyze the data, perform some plots by ourselves as well. Based on these analysis, we discuss about the result of the experiment and put forward some suggestions that may help the experiment to be better.

A MEASUREMENT UNCERTAINTY ANALYSIS

A.1 DEMONSTRATION OF MALUS' LAW

In this section, we record data as shown in Table 3. We want to calculate the uncertainty for I/I_0 and $\cos^2\theta$. The calculations will shown as follows:

- (1) The uncertainty of I/I_0 can be obtained through the following formula:

$$u_{I/I_0} = \sqrt{\left(\frac{I}{I_0^2}\right)^2 u_{I_0}^2 + \left(\frac{1}{I_0}\right)^2 u_I^2}$$

Take $\theta = 5^\circ$ as a sample calculation, we are able to obtain the uncertainty accordingly:

$$u_{I/I_0} = \sqrt{\left(\frac{1.203}{1.211^2}\right)^2 0.001^2 + \left(\frac{1}{1.211}\right)^2 0.001^2} = 0.0012$$

with relative uncertainty:

$$r_{u_{I/I_0}} = \frac{u_{I/I_0}}{I/I_0} \times 100\% = 0.12\%$$

Then, all the uncertainties are given in the following Table:

- (2) The uncertainty of $\cos^2\theta$ can be obtained through the following formula:

θ	I/I_0	u_{I/I_0}	$r_{u_{I/I_0}}$	θ	I/I_0	u_{I/I_0}	$r_{u_{I/I_0}}$
0°	1.0000	0.0012	0.12%	50°	0.4178	0.0009	0.2%
5°	0.9934	0.0012	0.12%	55°	0.3493	0.0009	0.3%
10°	0.9727	0.0012	0.12%	60°	0.2642	0.0009	0.3%
15°	0.9265	0.0011	0.12%	65°	0.1841	0.0008	0.4%
20°	0.9273	0.0011	0.12%	70°	0.1272	0.0008	0.6%
25°	0.8720	0.0011	0.13%	75°	0.0735	0.0008	1.1%
30°	0.7853	0.0010	0.13%	80°	0.0339	0.0008	2%
35°	0.7060	0.0010	0.14%	85°	0.0107	0.0008	7%
40°	0.6078	0.0009	0.15%	90°	0.0000	0.0008	/
45°	0.5120	0.0009	0.18%				

Table 11: Uncertainty of I/I_0 .

$$u_{\cos^2\theta} = \sqrt{\left(\frac{\partial \cos^2\theta}{\partial \theta}\right)^2 u_\theta^2} = 2\sin\theta \cos\theta u_\theta$$

Take $\theta = 5^\circ$ as a sample calculation, we are able to obtain the uncertainty accordingly:

$$u_{\cos^2\theta} = 2\sin 5^\circ \cos 5^\circ \frac{2}{180}\pi = 0.006$$

with relative uncertainty:

$$r_{u_{\cos^2\theta}} = \frac{u_{\cos^2\theta}}{\cos^2\theta} \times 100\% = 0.6\%$$

Then, all the uncertainties are given in the following Table:

θ	$\cos^2\theta$	$u_{\cos^2\theta}$	$r_{u_{\cos^2\theta}}$	θ	$\cos^2\theta$	$u_{\cos^2\theta}$	$r_{u_{\cos^2\theta}}$
0°	1	0	0	50°	0.41	0.03	7%
5°	0.992	0.006	0.6%	55°	0.33	0.03	9%
10°	0.970	0.012	1.2%	60°	0.25	0.03	12%
15°	0.933	0.017	1.8%	65°	0.18	0.03	17%
20°	0.88	0.02	2%	70°	0.12	0.02	17%
25°	0.82	0.03	4%	75°	0.067	0.017	25%
30°	0.75	0.03	4%	80°	0.030	0.012	40%
35°	0.67	0.03	4%	85°	0.008	0.006	75%
40°	0.59	0.03	5%	90°	0	0	/
45°	0.50	0.03	6%				

Table 12: Uncertainty of $\cos^2\theta$.

A.2 LINEARLY POLARIZED LIGHT AND THE HALF-WAVE PLATE

In this section, for the left column (Rotation angle of the 1/2-wave plate), since it is measured directly without calculations, what we need to consider is only the type-B uncertainty, i.e.

$$u_{\theta_1} = u_{\Delta_B} = 2^\circ$$

For the right column (Rotation angle of the analyzer), however, it is not the case, since we ought to calculate $\theta_2 = \Delta\theta$, i.e., we have the function:

$$\theta_2 = \theta_0 - \theta$$

Hence the uncertainty can be calculated as:

$$u_{\theta_2} = \sqrt{\left(\frac{\partial\theta_2}{\partial\theta_0}\right)^2 (u_{\theta_0})^2 + \left(\frac{\partial\theta_2}{\partial\theta}\right)^2 (u_\theta)^2} = \sqrt{(u_{\theta_0})^2 + (u_\theta)^2} = \sqrt{2^2 + 2^2} \approx 3^\circ$$

A.3 CIRCULARLY AND ELLIPTICALLY POLARIZED LIGHT AND THE 1/4-WAVE PLATE

Since $f(I, I_0) = \sqrt{I/I_0}$, we are able to find out its uncertainty through the following formula:

$$u_{\sqrt{\frac{I}{I_0}}} = \sqrt{\left(\frac{\partial f}{\partial I}\right)^2 (u_I)^2 + \left(\frac{\partial f}{\partial I_0}\right)^2 (u_{I_0})^2} = \sqrt{\frac{1}{4II_0} u_I^2 + \frac{I}{4I_0^3} u_{I_0}^2}$$

By the formula above, we are able to calculate all the uncertainties. Here we take the case when $I = 0.047 \mu A$ with $I_0 = 1.452 \mu A$ as an example. The uncertainties can be calculated as:

$$u_{\sqrt{\frac{I}{I_0}}} = \sqrt{\frac{1}{4 \times 0.047 \times 1.452} \times 0.001^2 + \frac{0.047}{4 \times 1.452^3} \times 0.001^2} = 0.002 \mu A$$

Now, all the uncertainties in this part can be calculated similarly.

A.3.1 ROTATION ANGLE 0°

Rotation angle of 1/4-wave plate: 0°					
Maximum Electric Current I_0			$1.452 \pm 0.001 \mu A$		
θ	$\sqrt{\frac{I}{I_0}}$	$u_{\sqrt{\frac{I}{I_0}}}$	θ	$\sqrt{\frac{I}{I_0}}$	$u_{\sqrt{\frac{I}{I_0}}}$
0°	0	/	180°	0.037	0.009
10°	0.180	0.002	190°	0.1874	0.0018
20°	0.3371	0.0010	200°	0.3452	0.0010
30°	0.5007	0.0007	210°	0.5163	0.0007
40°	0.6550	0.0006	220°	0.6243	0.0006
50°	0.7460	0.0005	230°	0.7339	0.0005
60°	0.8806	0.0005	240°	0.8724	0.0005
70°	0.9473	0.0005	250°	0.9488	0.0005
80°	0.9886	0.0005	260°	0.9872	0.0005
90°	1.0000	0.0005	270°	0.9997	0.0005
100°	0.9760	0.0005	280°	0.9721	0.0005
110°	0.9219	0.0005	290°	0.9170	0.0005
120°	0.8072	0.0005	300°	0.8093	0.0005
130°	0.7376	0.0005	310°	0.7404	0.0005
140°	0.5938	0.0006	320°	0.5909	0.0006
150°	0.4643	0.0008	330°	0.4621	0.0008
160°	0.3083	0.0011	340°	0.3049	0.0011
170°	0.1485	0.002	350°	0.166	0.002

Table 13: Measurement data for the 1/4-wave plate (Rotation angle 0°).

Based on our analysis above, we obtain the uncertainty of $\sqrt{I/I_0}$ when the angle is 0° .

A.3.2 ROTATION ANGLE 20°

Rotation angle of 1/4-wave plate: 20°					
Maximum Electric Current I_0			$1.201 \pm 0.001 \mu A$		
θ	$\sqrt{\frac{I}{I_0}}$	$u_{\sqrt{\frac{I}{I_0}}}$	θ	$\sqrt{\frac{I}{I_0}}$	$u_{\sqrt{\frac{I}{I_0}}}$
0°	0.3604	0.0011	180°	0.3707	0.0011
10°	0.4111	0.0010	190°	0.4414	0.0009
20°	0.5336	0.0008	200°	0.5467	0.0008
30°	0.6060	0.0007	210°	0.6342	0.0007
40°	0.7068	0.0007	220°	0.7306	0.0006
50°	0.7939	0.0006	230°	0.8121	0.0006
60°	0.8652	0.0006	240°	0.8767	0.0006
70°	0.9184	0.0006	250°	0.9355	0.0006
80°	0.9557	0.0006	260°	0.9743	0.0006
90°	0.9895	0.0006	270°	1.0000	0.0006
100°	0.9562	0.0006	280°	0.9794	0.0006
110°	0.8767	0.0006	290°	0.8992	0.0006
120°	0.7981	0.0006	300°	0.8268	0.0006
130°	0.7015	0.0006	310°	0.7185	0.0006
140°	0.5928	0.0007	320°	0.6202	0.0007
150°	0.4811	0.0009	330°	0.4905	0.0009
160°	0.4009	0.0011	340°	0.4386	0.0010
170°	0.3510	0.0012	350°	0.3673	0.0011

Table 14: Measurement data for the 1/4-wave plate (Rotation angle 20°).

Based on our analysis above, we obtain the uncertainty of $\sqrt{I/I_0}$ when the angle is 20°.

A.3.3 ROTATION ANGLE 45°

Rotation angle of 1/4-wave plate: 45°					
Maximum Electric Current I_0 0.705 ± 0.001 μA					
θ	$\sqrt{\frac{I}{I_0}}$	$u_{\sqrt{\frac{I}{I_0}}}$	θ	$\sqrt{\frac{I}{I_0}}$	$u_{\sqrt{\frac{I}{I_0}}}$
0°	0.9317	0.0010	180°	0.9438	0.0010
10°	0.9446	0.0010	190°	0.9400	0.0010
20°	0.9400	0.0010	200°	0.9416	0.0010
30°	0.9483	0.0010	210°	0.9506	0.0010
40°	0.9602	0.0010	220°	0.9617	0.0010
50°	0.9770	0.0010	230°	0.9785	0.0010
60°	0.9907	0.0010	240°	0.9893	0.0010
70°	0.9828	0.0010	250°	1.0000	0.0010
80°	0.9907	0.0010	260°	0.9957	0.0010
90°	0.9864	0.0010	270°	0.9950	0.0010
100°	0.9857	0.0010	280°	0.9843	0.0010
110°	0.9821	0.0010	290°	0.9763	0.0010
120°	0.9885	0.0010	300°	0.9734	0.0010
130°	0.9683	0.0010	310°	0.9668	0.0010
140°	0.9595	0.0010	320°	0.9572	0.0010
150°	0.9438	0.0010	330°	0.9438	0.0010
160°	0.9378	0.0010	340°	0.9355	0.0010
170°	0.9340	0.0010	350°	0.9378	0.0010

Table 15: Measurement data for the 1/4-wave plate (Rotation angle 45°).

Based on our analysis above, we obtain the uncertainty of $\sqrt{I/I_0}$ when the angle is 45°.

A.3.4 ROTATION ANGLE 70°

Since both the values are measured directly, the uncertainty is equal to their type-B uncertainty. Namely, we have:

$$u_{\theta} = 2^{\circ}$$

$$u_I = 0.001[\mu A]$$

B DATA SHEET

See the attached data sheet.

References

[1] Krzyzosiak, M. & VP141 TA Groups. *Exercise 5 - lab manual [rev 4.3].pdf*. 2019.